

NEW

EndoVue® 24"

Value Endoscopy Visualization



HD ENDOVue

- Optimized for endoscopy applications
- Outstanding image quality
- Multi-modality imaging capabilities
- Low power consumption and fanless cooling
- Robust and reliable
- Minimum cost of ownership

For the medical professional seeking to maximize their capital investment and upgrade their aging monitor, the EndoVue® series delivers the image performance of a high-end endoscopy display at a value price.

EndoVue is a versatile visualization platform, supporting multi-modality medical imaging requirements. It is plug-and-play compatible with endoscopic camera systems from leading manufacturers. Designed and built for medical use, it is FDA listed, TUV certified, and CE marked, offering patient safety, clinical performance, and durability compared to non-medical grade displays.

The EndoVue display accommodates standard analog and digital high-definition signals from a variety of medical imaging sources, including ultrasound, PACS, and patient monitor vital signs. It features a fanless cooling design to minimize the risk of spreading airborne contaminants within the sterile field, and industry-leading LED backlight technology for reduced power consumption.

EndoVue® 24"

Value Endoscopy Visualization



SPECIFICATIONS

| | |
|-----------------------|---|
| Screen Size | 24" diagonal |
| Resolution | 1920 x 1200 |
| Aspect Ratio | 16:10 |
| View Angle | 178° |
| Luminance | 300 cd/m ² typical |
| Color Gamut | BT.709 or SMPTE-C |
| Display Warranty | 1 Year |
| Dimensions (WxHxD) | 21.9" x 16.1" x 2.8"/566 x 408 x 70 mm |
| Weight | 14.7 lbs/6.7 kg |
| Green Compliance | REACH, RoHS-2, WEEE |
| Regulatory Compliance | ANSI/AAMI ES60601-1, CAN/CSA C22.2 No. 60601-1, CAN ICES-3(B)/NMB-3(B), FCC Class B, EN60601-1, EN60601-1-2, CE, MDD 93/42/EEC, CCC, EAC, Conflict Minerals, Class I Medical Device |

ORDER INFORMATION

EndoVue 24

| | |
|---------|---------------------------------------|
| 90K0070 | EndoVue 24" - North America |
| 90K0071 | EndoVue 24" - South America |
| 90K0072 | EndoVue 24" - Japan |
| 90K0073 | EndoVue 24" - China |
| 90K0074 | EndoVue 24" - Emerging Market |
| 90K0075 | EndoVue 24" - Europe |
| 26B0102 | EndoVue Display Stand - 19", 21", 24" |

3061 NW 82ND Ave • Miami, FL 33122
Ph.: 305-406-2222 • Fax: 305-406-2213
www.shmedical.com