

# LED-2000 50W Solid State LED Light Source



*Premium LED-2000 with 4-port turret*

## LED-2000 Medical Grade 50W LED Light Source

- Compact, efficient, durable sources of brilliant “daylight” quality light.
- Premium LED-2000 utilizes new high efficiency design, with no change in form, generating 50% greater outputs for more demanding applications.
- Color temperature: 6500°K nominal
- Typical LED life: >50,000 hours
- Low power consumption: 50W typical.
- Electronic (PWM) intensity control from 0-100%.
- Multi-port turret and single port versions available.
- Standard ports are ACMI, Olympus, Storz and Wolf; other ports available.
- Accepts 100-240V AC, 50/60 Hz input power for worldwide operation.
- Regulatory: UL60601-1, CAN/CSA C22.2 No. 601.1, EN60601-1-2; marked.
- Dimensions: 4.65” W x 3.32” H X 7.25” D (~10” overall)
- Weight : 2.5 lbs (1.14 kg)