

NEW



## EndoVue® 32"

### Value Endoscopy Visualization



- Optimized for endoscopy applications
- Outstanding image quality
- Multi-modality imaging capabilities
- Low power consumption and fanless cooling
- Robust and reliable
- Minimum cost of ownership

For the medical professional seeking to maximize their capital investment and upgrade their aging monitor, the EndoVue® series delivers the image performance of a high-end endoscopy display at a value price.

EndoVue is a versatile visualization platform, supporting multi-modality medical imaging requirements. It is plug-and-play compatible with endoscopic camera systems from leading manufacturers. Designed and built for medical use, it is FDA listed, TUV certified, and CE marked, offering patient safety, clinical performance, and durability compared to non-medical grade displays.

The EndoVue display accommodates standard analog and digital high-definition signals from a variety of medical imaging sources, including ultrasound, PACS, and patient monitor vital signs. It features a fanless cooling design to minimize the risk of spreading airborne contaminants within the sterile field, and industry-leading LED backlight technology for reduced power consumption.

# EndoVue® 32"

## Value Endoscopy Visualization



### SPECIFICATIONS

Screen Size	31.5" diagonal
Resolution	1920 x 1080
Aspect Ratio	16:9
View Angle	178°
Luminance	400 cd/m <sup>2</sup> typical
Color Gamut	BT.709 or SMPTE-C
Display Warranty	1 Year
Dimensions (WxHxD)	29.7" x 19.1" x 3.9"/753 x 485 x 98 mm
Weight	21.3 lbs/9.7 kg
Green Compliance	REACH, RoHS-2, WEEE
Regulatory Compliance	ANSI/AAMI ES60601-1, CAN/CSA C22.2 No. 60601-1, CAN ICES-3(B)/NMB-3(B), FCC Class B, EN60601-1, EN60601-1-2, CE, MDD 93/42/EEC, CCC, Conflict Minerals, Class I Medical Device

### ORDER INFORMATION

#### EndoVue 32

90K0080	EndoVue 32" - North America
90K0081	EndoVue 32" - South America
90K0082	EndoVue 32" - Japan
90K0083	EndoVue 32" - China
90K0084	EndoVue 32" - Emerging Market
90K0085	EndoVue 32" - Europe



**3061 NW 82<sup>ND</sup> Ave • Miami, FL 33122**  
**Ph.: 305-406-2222 • Fax: 305-406-2213**  
**[www.shmedical.com](http://www.shmedical.com)**